















USER MANUAL

DISCO/ DISCO E-BIKE/ DISCO PEDELEC













TABLE OF CONTENTS

1	Foreword:	4
2	Introductory information:	6
2.1	Use	6
2.2	Accessibility	6
2.3	Adjustments	6
2.4	Equipment	6
2.5	Transportation  	6
2.6	Cleaning and maintenance	7
2.7	Auxiliary electric motor  	7
2.8	Technical specifications	7
3	Before using the bicycle for the first time:	8
3.1	Statutory provisions	8
3.2	Preparing the bicycle	8
4	Riding your PF bicycle	9
4.1	Riding instructions	9
4.2	Before taking your bicycle out for a ride	10
5	Safety instructions	10
5.1	Wear and tear of wheels and tyres	10
5.2	Warnings	10
6	Settings and adjustment:	11
6.1	Handlebars:	11
6.2	Seat:	11
6.3	Brakes	12
6.4	Gear	12
6.5	Adjustment of chains	13
7	Electrical system and motor  	14
7.1	Electric bicycle, Pedelec or E-bike?	14
7.2	Safety instruction  	14
7.3	Hazards for special user groups	15
7.4	Battery  	15
7.5	Charging the battery  	16
7.6	Using the electric bicycle	17
8	Accessories	21



9	Maintaining your bicycle	22
9.1	Important  	22
9.2	Daily check	23
9.3	Periodic maintenance – and first service check!	24
9.4	Cleaning  	25
9.5	Lubrication	25
9.6	Tightening of bolts, etc.	25
9.7	Repair and replacement	26
10	Disposal  	26
11	Warranty and obligations	27
12	Error codes and troubleshooting  	27
13	EC Declaration of Conformity	31



1 FOREWORD:

Congratulations on your new PF bicycle. We hope that you will enjoy riding it for many years to come! Please read carefully this User and Operating Manual. To be able to take maximum advantage of your bicycle, it is important that you adjust it properly and that you follow the advice we have for you in this User and Operating Manual. Take particular note of Section 4.



Riding your PF bicycle

IMPORTANT: For safety reasons, it is important that the bicycle be serviced at all times by a person with professional training and insight into bicycle technology and function. Take particular note of Section 9 *Maintaining your bicycle*.

PF mobility has an extensive range of accessories that can contribute to improving the use and experience of your PF bicycle. Our accessories are listed on our website, but a selection of them is also presented under Section 8 *Accessories*.


If you have any further questions about the use of your PF bicycle or this Manual, you are more than welcome to contact us at:

PF mobility ApS
Bjerregårdvej 15, Timring
DK – 7480 Vildbjerg
Tel.: +45 99 92 06 00
Fax: +45 99 92 06 05

Mail: info@pfmobility.dk
www.pfmobility.dk

With reservations for misprints and changes to this User Manual.

Product details:	
HMI - Number	
Frame number	
Lock number	
Controller keys	
Battery number	
Delivered on	



Fill in this table as per the specifications in your invoice.
The frame number is also the bicycle's unique ID number. The frame number is placed on the side of the stem – see image.



2 INTRODUCTORY INFORMATION:

2.1 Use


Your PF bicycle has been designed for people with functional disabilities – whether physical or mental.

These include problems with balance, in the muscular system, etc.

The bicycle is designed for use on the common public roads and bicycle paths. Always ride the bicycle on an even solid surface and adapt your speed to the conditions at hand. Comply with the provisions of the Danish Road Traffic Act. PF mobility may not be held liable for accidents resulting from improper use or for damage resulting from lack of maintenance.

2.2 Accessibility

We recommend that people with cognitive functional disabilities, visual impairments or reading disabilities, who cannot read and understand this Manual, turn to the dealer where they bought their PF bicycle. The dealer will review the manual and reply to any questions they may have.



2.3 Adjustments

To benefit fully from your new PF bicycle, we ask you to carefully peruse this Manual. Take note of the sections regarding settings and adjustments. Your dealer will help you with any questions about specific settings. Even if the bicycle is made of quality components, it still needs maintenance. You can easily make many of the minor adjustments yourself. This will both extend the service life of the bicycle and make the bicycle ride true enjoyment.

2.4 Equipment

There is a large range of accessories designed for your bicycle. See Section 8 regarding accessories or visit our website.

2.5 Transportation

Fasten the bicycle securely if you need to transport it in a car or in a trailer. When you use hose clamps and similar, make sure that they do not jam any chains, wires and cables. During transportation, your bicycle may be exposed to moisture, dirt and other aggressive substances. This can affect the electrical system, which is why you must cover the bicycle with plastic or tarpaulin when you transport it in an open trailer.

During transportation in a car or trailer:



- Remove the battery and store it in a cool and dry place.
- In humid weather, protect the bicycle by covering it with plastic or tarpaulin.



Note: Do NOT ship the Li-Ion battery as an ordinary package. Li-Ion batteries are considered hazardous goods and should be marked in accordance with the instructions in the ADR convention (UN3840).





2.6 Cleaning and maintenance

Your bicycle needs day-to-day maintenance.

We recommend that you have an annual maintenance check performed by a bicycle professional or PF mobility. A review of maintenance checkpoints and methods of cleaning is available under Section 9 *Maintaining your bicycle*.






2.7 Auxiliary electric motor

For bicycles equipped with an electric motor, kindly refer to the *E-Bike and Pedelec* section.

Sections containing important information for you as an owner of an e-bike are marked in the heading with  

These sections contain warnings you should have a special focus on. Make sure you abide by these.

The following symbols are used in this Manual:

 	Warning: This symbol warns about a hazard to your health and points at a potential risk of injuries.
	Note: Indicates possible material damage
	Important: Contains general instructions for safe use and special technical features or directions.
	Information: This symbol refers to tips or special information.

2.8 Technical specifications

	Model	DISCO	DISCO MEDI	DISCO SMALL
Data				
Lenght cm		185	171	154
Width cm		75	68	68
Weight (w/o accessories)		Ca. 44	Ca. 29	Ca. 28
Extra weight E-Bike/ Pedelec		Ca. 11	Ca. 11	Ca. 11
Maximum user weight kg		125	125	125 kg
Maximum load kg		150	150	150 kg
Inside leg (Min)		55	45	40 cm
Mounting/Dismounting height		26	25	20 cm
Wheel size		20	16	16
Max speed motor km/t E-bike / Pedelec		15/ 25	15/ 25	15/ 25



EPAC and Pedelec models

Steering		
Operating voltage		36V
Max amperage		Depending on version and use, 10 to 30A
Motor		
Nominal voltage		22.8 V
Motor output		250W (DIN EN 60034-1)
Nominal/Peak momentum		11.8 Nm / up to 60 Nm
Battery		
Type		Lithium-Ion battery
Nominal voltage		36 V
Capacity		11 Ah, 400 Wh (Rear Rack – standard) 14 Ah, 500 Wh (Downtube – option)
Charging period		Approx. 6 hours at charging current of 2 A (Rear Rack) Approx. 7 hours at charging current of 2 A (Downtube)
Temperature limitations	Normal operation	-10 to 45 °C
	Charging	10 to 35 °C
	Storage	-10 to 45 °C

3 BEFORE USING THE BICYCLE FOR THE FIRST TIME:

3.1 Statutory provisions

The bicycle is equipped from the factory with statutory equipment pursuant to the Danish Executive Order regarding Design and Equipment, etc. of Bicycles:

We recommend that the bicycle is always fitted with lights as per the applicable legislation. Consult PF mobility or your dealer.

Even if riding with a bicycle helmet is not compulsory by law, we strongly recommend that you always wear a helmet.

It is your responsibility to make sure that defective parts are replaced so that the bicycle remains in lawful condition and thereby safe to ride at all times. PF mobility recommends that you only use original spare parts that ensure that the bicycle meets the applicable legislation.

3.2 Preparing the bicycle

The bicycle is delivered assembled from the factory. If the bicycle has been received from a carrier/forwarding agent, it may be necessary to adjust and fasten properly the handlebars and the seat. See Section 2.3 *Adjustments*.

If your bicycle has been delivered by a consultant from PF mobility, they will review all necessary adjustments together with you. You are subsequently always welcome to contact your dealer or PF mobility for any additional help or instructions.



4 RIDING YOUR PF BICYCLE

To get as much riding joy and benefit as much as possible from your PF bicycle for many years, it is important that you use the bicycle properly.



You can ride the bicycle on different surfaces. PF mobility recommends that you ride on even solid surfaces (preferably asphalt or slabs).

The bicycle is not designed for riding in rough terrain or on a soft surface as this can overload the structure and the frame or the electrical system.

4.1 Riding instructions

Avoid:

- engaging the brakes at a high speed as this involves a risk of skidding and overturning;
- passing over holes and tracks as this can cause you to lose your balance and overturn with the bicycle;
- forcing kerbs – use the ramps instead. If this, nevertheless, is necessary, get off the bicycle and pull the bicycle up or down over the edge. If you have a bicycle with Direct Power, you can use the auxiliary motor's help function. See section 7.6.4 *Start and help function*.

Adapt your speed:

- when riding in rough terrain. The bicycle will behave differently in such conditions, which is why you should adjust your speed to avoid unpleasant and outright dangerous situations;
- in bends. If you go around a bend at too high a speed, you risk overturning;
- to the weather and other conditions. Riding on wet and slippery roads poses a particularly high risk, which is why you should ride carefully and at a reduced speed – in particular in curves and bends as well as when riding downhill.

Note that:

- the bicycle is wide – watch out for kerbstones, roadsides, chicanes, people and other obstacles;
- the bicycle can make undesired movements while you mount or dismount it. ALWAYS use the parking brake when you park the bicycle and, in particular, when you mount or dismount the bicycle.

If you are also carrying luggage / wheelchair / walker, make sure that it is securely fastened. You can find the right accessory to ensure secure transportation and fastening of your bicycle in the section about accessories.



4.2 Before taking your bicycle out for a ride

Check daily, before each ride:

- ✓ That the bicycle is in a perfect condition;
- ✓ The tyre pressure;
- ✓ That the light and statutory reflectors are installed.
- ✓ The following brakes: Handbrake in front, handbrake at the back as well as foot brake (assuming the model has a foot brake);
- ✓ Parking brake
- ✓ That the handlebars are properly fastened;
- ✓ That the seat is properly fastened;
- ✓ That there is no gap in the wheels;
- ✓ That the chains are lubricated and adjusted;
- ✓ That the bicycle is fitted with the statutory reflectors;
- ✓ That the bicycle is fully charged (if it is an E-Bike / Pedelec).



5 SAFETY INSTRUCTIONS

5.1 Wear and tear of wheels and tyres

You should perform regular (weekly) checks of the bicycle's tyres. If you see any hair cracks starting to form on the side of the tyre, replace the tyre. New tyres can be ordered from PF mobility or your dealer. We recommend that you always use puncture-proof tyres. The bicycle is delivered by default with tyres with high puncture protection and reflective edges.

Also check the wheel rims for cracks or ruptures. Bicycles that ride on very rough / uneven terrain are more receptive to wear and tear.



Defective wheel rims can lead to accidents and should therefore be replaced immediately!!!

5.2 Warnings

Take off any loosely hanging clothes and put on something that fits snugly and cannot get caught in rotating parts such as wheels and chains before you mount your PF bicycle.

Do not use your PF bicycle if any of its safety guards, e.g. wheel guards, chain guards or the guard over the rear axle are defective or missing.

Take extra care when adjusting the seat or armrests as inattention may cause you to jam your fingers. Only make adjustments when the bicycle is completely still.



6 SETTINGS AND ADJUSTMENT:

The following pages contain instructions for the daily adjustment options of the PF bicycle. Wherever you need to use tools, it is important after the adjustment to tighten all bolts, etc. you have loosened.

6.1 Handlebars:

Adjust the angle of the handlebars by loosening bolts "1" and then place the handlebars in the desired position. Tighten up bolts "1" with a torque of 10 N.m.

Adjust the height of the handlebars by loosening bolt "2" until the little plate lets loose of the cogs. Then you can change the height. Tighten back the bolt until the plate locks into the cogs.



6.2 Seat:

You can regulate the distance from the seat to the pedals by lifting release arm "1" underneath the seat frame. You will need to loosen the clamping lever "2" if you need to dismantle or adjust the seat vertically. When raising the seat height, you MUST fit the original spacer pipe to ensure that the seat remains at the required height (optional equipment).

There is a screw on the opposite side of the release arm that makes sure that the seat cannot slide off the rails. Check regularly that this screw is tight.



Do not raise the seat above the marking.

- use special spacer / extension pipe for the seat if you need to raise the seat (optional equipment – see list).



6.3 Brakes

All of our bicycles are fitted with at least 2 brakes.

The left handbrake handle operates the V-brake on the front wheel, whereas the right handbrake handle operates the brake disc fitted to the differential gear on the rear side of the bicycle. The bicycle can also be fitted with a foot brake (gear hub) or a freewheel wherever a foot brake is not an option.

The parking brake is fitted to the handlebars. Turn it up to activate the brake. To disengage the brake, give the handle a gentle pull until the lever falls down the low mark.

The brakes are in need of continuous adjustments. Turn the long adjusting screw marked in the picture counter-clockwise and fix the check nut farthest in towards the brake handle. Extensive adjustments of the V-brake should only be conducted by a professional.



6.4 Gear

Your PF bicycle is delivered by default with 7 gears made by Shimano.

Twist the rotating handle "1" with the thumb and index finger to switch gear.

Keep the pedals still when you change gear while riding.

It is necessary to make subsequent adjustments of the gear. You can do this, as follows: set the rotating handle "1" to gear 4. Then adjust the gear with the adjusting screw "2" until the yellow dots "3" on the gear hub are parallel.




6.5 Adjustment of chains

Front chain (A): Loosen the hex bolt on the left side with the 8-mm hex key (1). Turn the double gearwheel clockwise with the 17-mm combination wrench (1) to tighten the chain. Hold the bolt of the double gearwheel in position with the 17-mm combination wrench (1) to adjust the chain and then tighten the bolt securely to (1).

Middle chain (B) and back chain (C): Loosen the hub nuts and the brake lever's nut and screw. Loosen the counter nut of the adjusting screw using the 17-mm fork spanner (3). You can now adjust the medium and back chain with the adjusting screw and the chain adjusters.

Remember to tighten well after each adjustment!!!





7 ELECTRICAL SYSTEM AND MOTOR

Peruse this chapter carefully if your bicycle is equipped with an auxiliary electric motor. Save these instructions and remember to bring them if you take the bicycles to other repairmen.

Non-compliance with these instructions can lead to personal injuries or damage to the bicycle or its components.

PF mobility may not be held liable for any damage or loss resulting from failure to comply with these instructions.

7.1 Electric bicycle, Pedelec or E-bike?

Electric bicycle is a common designation that is used for both a Pedelec and an E-bike. The auxiliary electric motor of a Pedelec is only activated when you are pushing the pedals. The bicycle's auxiliary motor is only active at speeds of up to 25 km per hour. Nevertheless, the bicycle can be set to cut out the auxiliary motor at a lower speed – see your order confirmation. The programmed "max" speed is also stated on the programming label of the bicycle that you will find on the steering behind the battery.



An E-bike's auxiliary motor can be activating even if you are not pushing the pedals. The bicycle is fitted with a throttle that determines the motor's output – however, only up to a speed of max. 15 km per hour, when the motor is cuts out. Things can be different – see the same label as above.

Both models have a starting aid of up to 6 km per hour.

7.2 Safety instruction



Risk of breakdown of load-bearing parts.

A defective motor can cause damage to load-bearing parts and bearings. In turn, this can cause you to crash or fall.

- Stop using the bicycle at once and contact without delay PF mobility or a bicycle expert in order to replace the defective motor.



Hazard of fault or fire

Defective electrical connections or cables can cause short circuits. This can cause electric faults or, in the worst case scenario, fire!

- Replace any ruined or defective electrical connections and components.



Hazard of unintended operation of the motor and other rotating parts

The motor can be activated inadvertently if the electrical system remains "on" while your bicycle is undergoing maintenance. Clothes can get caught into the rotating parts and cause personal injuries.

- Always remove the battery from the bicycle before you start any work on the bicycle.

7.3 Hazards for special user groups

Children and adolescents are allowed to use electric bicycles on the public road network, but only under the supervision of an adult.

The bicycle should be parked out of the reach of children or persons who are not capable of handling the specified risks.

Riding an electric bicycle is different than riding an ordinary bicycle. This is why we recommend that you practise in a place without traffic before you ride in traffic.

7.4 Battery

A Li-Ion battery combines a low weight with a very high current capacity. This makes it rather compact and contributes to keeping the weight of the bicycle low.

Your PF bicycle is delivered by default with a "Rear Rack" battery, but it can also be delivered with a "Downtube" battery (optional equipment) that has extra battery capacity.



"Rear Rack"



"Downtube"

Li-Ion batteries may only be charged with a special charging system. You can significantly increase the service life of the battery by ensuring proper recharging and preventing the battery from discharging completely. To prevent this and help you as owner of the bicycle, we have taken these factors into account in the enclosed charger. For the same reason, you may only use the charger to charge the battery of your bicycle.



Fire hazard

Incorrect recharging of the Li-Ion battery may cause the battery or charger to heat substantially. This can cause a fire!

- Always use the accompanying charger when recharging the battery. The charger is NOT suited for outdoor use.
- Before you connect your charger to the mains, you must ensure that the mains voltage corresponds to the permitted voltage of the charger. The voltage is stated on the rating plate of the charger.
- Only charge the Li-Ion battery in a dry, non-flammable environment.

Mechanical damage to the Li-Ion battery or the charger can cause malfunction and short-circuit. In turn, this can cause fire.

- Any form of manipulation of the battery or charger is strictly forbidden.

- **Immediately** replace the battery if ruined and dispose of it safely, e.g. by handing it over to a municipal receiving station.



Risk of electric shock

Chargers with defective cords or connectors can cause an electric shock.

- Never plug defective connectors or cables to the mains.
- Immediately replace any defective components and cables.
- Protect the charger against penetration of water and moisture. If water has penetrated into the charger, pull out the connector and have the charger inspected by a professional.
- Condensation can form if the charger has been exposed to a sudden temperature increase. If this happens, you must wait until the charger has the same temperature as its surroundings before you connect the charger to the mains. Always store the charger wherever it is used.
- Only use the charger to charge the supplied Li-Ion battery. Other use is not permitted.

7.5 Charging the battery

Be sure to fully charge the battery before you use your bicycle for the first time. The battery can be charged while it is installed to the bicycle. The battery can also be detached and charged in a more appropriate location.

1. Open the dust cap and insert the charging plug.



Støvhætte – Downtube

Støvhætte – Rear rack

To ensure proper connection, the connection between the charger and the battery is designed with a groove and a protrusion. These are marked with an arrow on both the battery and the charger. Make sure that you have positioned the charging plug correctly – do not apply force.

2. Connect the charger to a power socket.
3. Once the battery is charged, the red indicator on the charger will light and the green indicators on the battery will flash (only Rear Rack).
4. When the charger's indicator lights permanent green, the battery will be charged, and you can disconnect the power.



The charging time is approx. 6 hours for full charging of the Rear Rack battery and approx. 7 hours for the Downtube battery. The charger then switches to maintenance charging, which ensures that the battery does not discharge. If the battery is not fully charged, the bicycle's range will be reduced. Do not charge the battery at temperatures below 10 °C or over 35°C, in direct sunlight, next to radiators and similar. When the temperature outside is below 0°C, we recommend that you charge your battery in a heated room. If you do not use the bicycle for long periods of time, you should remove and store the battery in a dry, frost-proof place. Make sure to fully charge the battery before storage and again before you are about to use the bicycle.



Every time you install the battery on the bicycle, make sure that the casing is locked.

7.6 Using the electric bicycle

7.6.1 Instructions for use




Risk of hot surfaces

The battery and the motor can heat after long and demanding rides.

- Do not touch the battery or the motor after very demanding rides.

7.6.2 Operation of the electric bicycle

Your electric bicycle is equipped with a displayed and an operating switch.



- A. Assist / output level – determines the performance of the electric motor.
- B. Bicycle light indicator – lights when the light of the bicycle is on (optional equipment).
- C. Battery level – how much power there is left.
- D. Speedometer showing your current speed – as well as error codes (see section 12).
- E. Function display for different functions (see Section 7.6.7)
- F. Display line – displays the value of the individual functions selected in "E" as well as "Error" (see Section 12).



There is an operating switch with three buttons on the handlebars: MODE, "Up" and "Down" arrow.
All functions are operated via these buttons.

7.6.3 Turning the bicycle on and off



Press and hold the "MODE" button for three seconds until the display lights.
To turn off the system, press and hold the "MODE" button for 2 seconds.



The bicycle's system turns off automatically when the bicycle has stood still for 10 minutes.

WARNING:



Do **not** place your feet on the pedals when you turn on the bicycle as you risk to inadvertently set the bicycle in motion or to generate a system error.

- Mount the bicycle, but keep your feet on the ground.
- Turn on the bicycle's electrical system.
- Place your feet on the pedals and ride.

7.6.4 Start and help function

Start is a function that activates the motor in order to help you get going if you cannot do it yourself with the pedals. This function will help you achieve a speed of up to 6 km/h. The function can also be used for help with towage when you are pushing the bicycle yourself.






To activate the start help, push the "Up" button for approx. 1 second and then hold it in for as long as it is necessary.
The start help will cut out once you release the button.

For models with a throttle, e.g. an E-bike, the start help is activated by turning the throttle.

7.6.5 Adjustment of Assist level

The electrical system can be adjusted and set to three different assist / output levels. The higher output level you opt for, the greater the motor's output. The selected and relevant output level is visible in the top right corner of the display.

Assist / output level			
Motor output	~ 30-35 %	~ 65%	100 %
Start help	•	••	•••
Range	•••	••	•



To select a higher assist / output level, give the "Up" button one short push.



To select a lower assist / output level, give the "Down" button one short push.



Irrespective of which assist / output level you choose, the auxiliary motor's pre-programmed max speed will remain the same – See Section 7.1 *Electric bicycle, Pedelec or E-bike?*.

7.6.6 Power level

You can always see how much power there is in the battery in the top right corner of the display under "Battery". The display indicates the power level with 0 to 6 bars where 6 bars (solid coloured) indicates that the battery is fully charged and a single flashing bar indicates that the battery is almost dead.



If there are no bars in the display, the battery is empty, and the system will shut down before long. After that you will not be able to use any of the electric functions of the bicycle.

It will first be possible to use the electric functions of the bicycle once the battery is charged again.

If you charge the battery while it is installed on your bicycle, you will also be able to read off the charging status of the battery in the display.



WARNING:

Make sure that the battery is charged enough so that you have enough juice for the return trip.

7.6.7 Functions

The display has different display functions you can use while you are riding your bicycle. These are:

- Dist: Displays the overall number of km you have ridden.



- Trip: Trip counter displaying the number of km ridden since last reset.
- Time: Journey time:
- AVG: Average speed.
- EstT: Displays how long you can expect to ride at the current battery level, use and terrain – the value is only indicative!
- EstD: Displays the distance you can expect to cover at the current battery level, use and terrain – the value is only indicative!
- PIN: Not active for this model.



You can switch among individual display functions by pressing the "MODE" button several times. Each push switches over to the next display function.

7.6.8 Light

Your electric bicycle can be equipped with a light operated via the bicycle's electrical system (optional equipment).



If you have purchased an optional light connected to the electric bicycle's electrical system, you can turn it on by pressing the "Down" button for 2 seconds.

To turn off the light, press the "Down" button again for 2 seconds.



When the light is on, you will see a "Light" symbol at the top of the display.



8 ACCESSORIES

You can buy additional accessories and equipment for your PF bicycle. Get in touch with your dealer or PF mobility if you want to know more about your options.

Foot supports:



Footplate 5 sizes.



Pedal with toe clip



Pedal with toe clip and side support



Footplate with calf support



Training pedal



Pedal extender

Saddles & belts:



Armrest



Swivel seat



Junior seat



Seat console for junior seat



H-harness with holder



Extension pipe



H-harness with holder and head support



Harness



Back support with harness



Abduction block



Seat cover

Other accessories:



Basket



Topbox



Holder for walker



Holder for walking stick



Bike cover




Safety vest



Battery light



Dynamo light set



Bike computer

Lights:

Fits for:

- Disco
- Stabilo

- Duo
- PF Duo

9 MAINTAINING YOUR BICYCLE

9.1 Important

WARNING – Risk of malfunction in the event of incorrect or missing maintenance.



Incorrect or missing maintenance of an electrically operated bicycle can cause faults and damage to electrical and mechanical components! In this worst case scenario, this can cause you to crash and suffer a personal injury.

- The electrical system may only be maintained by a professional.



Hazard of unintended operation of the motor and other rotating parts

The motor can be activated inadvertently if the electrical system remains "on" while your bicycle is undergoing maintenance. Clothes, etc. can get caught into the rotating parts and cause personal injuries.

- Always remove the battery from the bicycle before you start any work on the bicycle.

Provided that they are both used as described and intended in this Manual, the electric motor and the battery are maintenance-free.

Have a professional conduct a safety check of the following items every 1,000 km:

- Intact and secure fastening of cables and components
- Review and testing of the complete electrical system.
- Safe use of the battery

9.2 Daily check

As mentioned above in Section 4.2 *Before taking your bicycle out for a ride*, you should check:

- Brakes and parking brake
- Chains
- Tyres and wheel rims
- Light and reflectors
- Battery condition

Your bicycle is equipped with a **parking brake** that is designed to prevent the bicycle from making any undesired movements during parking. This is particularly important when you mount or dismount the bicycle. You should test the parking brake and make sure that it operates adequately every day. We recommend that you always use the parking brake when you get off the bicycle. Adjust the brake when it does not operate adequately or as intended. For adjustments, see Section 2.3 *Adjustments*.



Test the **foot brake and handbrakes** daily to ensure that you can bring the bicycle to a halt quickly and safely.

Note that you can experience longer braking distance in humid weather.

Chain

To be able to use your bicycle safely and without any risk, you must check that the chain is properly adjusted and lubricated.

For the proper adjustment of brakes or chains, refer to a professional – see also Section 2.3 *Adjustments*.

Tyres and wheel rims



We recommend that you check daily – or at least weekly – that the tyres are properly inflated and not worn. The correct tyre pressure is 4-5 bar.



Since the bicycle uses rim brakes, you should also check the wheel rim of the front wheel for wear and tear and cracks. For more information, refer to Section 5.1.

Light and reflectors



Check that the bicycle's reflectors are whole and clean and that all statutory reflectors are fitted to the bicycle. White in front, red at the back, yellow on the pedals and white again on the tyres (standard equipment) or yellow on the spokes.

If you ride the bicycle at night, before you mount, you must check that the **light** of your bicycle works.

  Before you leave home, you must also check that the **battery** is charged – see Section 7.5

9.3 Periodic maintenance – and first service check!



In order to be in a good and safe mechanical condition, your bicycle is in need of periodic maintenance. The maintenance inspection should cover, as a minimum:

- Brakes
- Tyres and wheel rims
- All moving parts and bearings, crank, pedals, etc. (clearance and mobility)
- Chain
- Gear
- Electrical system and battery
- Frame
- Tightening of bolts, screws and nuts

The inspection should be conducted by a professional. We recommend a maintenance inspection every 1,000 km or at least once a year. You can buy an inspection from PF mobility through your sales consultant.

Service check



A service check must be carried out two months after the purchase of each new bicycle, cf. Section 9.3.

The gear, wires, bolts, chains, spokes, etc. of each new bicycle need some time to "settle", which is why it is necessary to check carefully and tighten / adjust individual parts.

This inspection must be conducted by a professional!

Lack of maintenance / service inspection can cause hazardous incidents and void the warranty.

9.4 Cleaning

Important information for electric bicycles.



- Remove the battery before you start cleaning your bicycle.

You can clean the bicycle with a moist – but not wet – soft cloth and water with a bit of detergent added, e.g. bicycle shampoo or similar.



Never use high-pressure cleaner or water jet.

Do not use a hot water or steam cleaner, a high-pressure cleaner or strong water jet / spray to clean your bicycle. If water penetrates inside, it can ruin the electrical components, motor bearings and other moving parts of your bicycle.



9.5 Lubrication



Clean and lubricate the chain of your bicycle, as needed. The more often you use the bicycle, the more often the chain will be in need of lubrication. Riding in humid and wet weather increases the need for lubrication. We recommend that the chain be lubricated at least once a month and more often in the winter season. You can buy a suitable detergent and lubricant from your bicycle professional. As a rule, none of the other moving parts is in need of regular lubrication, but should only be inspected periodically by a professional and lubricated, as needed.

9.6 Tightening of bolts, etc.



It is important to ensure that all critical components of the bicycle are attached properly. These parts need to be tightened with the proper torque. The required tightening force is stated in the following table:



Component:	Torque (+/- 5 Nm)
Wheel nut, front	45 Nm
Wheel nut, back	45 Nm
Stem – height	Indicated on the stem
Exocentric disc (front sprocket)	70 Nm
Rear beam	35 Nm
Rear axle – bearing	20 Nm
Saddle	40 Nm

9.7 Repair and replacement

To guarantee your safety, PF mobility recommends that you only use original spare parts. This will also allow you to benefit the most from the use of the bicycle. PF mobility has spare parts on stock and can deliver parts from one day to another. Contact your dealer or PF mobility for assistance and delivery of parts as well as replacement, if relevant.

10 DISPOSAL



Do not dispose of electronic and electrical components as ordinary household waste as this is prohibited by law. You as a user should therefore make sure that all electronic and electrical components are disposed of via an approved public waste management scheme, through a dealer or by PF mobility.

By complying with this and recycling or using an approved return scheme, you will contribute to making a real difference and help protect the environment.



Dispose of electrical components in accordance with the applicable national legislation.

For advice about disposal, you can contact the local authorities or your dealer.



11 WARRANTY AND OBLIGATIONS

All of PF's bicycles sold to individuals come with 2 years of warranty; cf. the provisions of the Danish Sale of Goods Act. In addition, there is 5-year warranty for the frame provided that the product has been used as described in this Manual.

Assuming ordinary use and periodic service, the expected service life of this bicycle is 5 to 10 years.


PF mobility may not be held liable for any damage to the bicycle itself or other objects or for injuries to people:

- caused by the abnormal use of the bicycle, **e.g.** in races or for riding in rugged terrain, etc.;
- if the user has disregarded or failed to comply with the instructions in this Manual;
- caused by inadvertent use of the PF bicycle;
- if there are any alterations to the structure of the bicycle without the producer's written consent;
- if the bicycle has been exposed to abnormal wear and tear or an overload;
- if the bicycle is not maintained normally or is exposed to extreme corrosion.

12 ERROR CODES AND TROUBLESHOOTING

In case there are problems or faults in the electrical system of your Direct Power electric bicycle, the display will show "E r r o r" along with an error number.

The different error codes are listed in the table below together with a description of possible causes and an appropriate solution proposal.

Error number	Error description and solution proposal
	Most errors can be solved by rebooting the system: Turn off the system, wait a bit and then turn it on again. If the error does not disappear after you have tried the proposed solution or occurs periodically, you need to contact your dealer or PF mobility.
1 <i>Error</i> Solution	Error Hardware Brake: <i>Cut out because of overvoltage</i> Remove the battery and check the connections. If the connections are in order, insert the battery again and reboot the electrical system (turn off – wait a bit – turn on again). Contact your dealer if the battery or the connections are damaged.
2 <i>Error</i> Solution	Error Hall <i>Signal between motor and steering interrupted.</i> Check the motor cable for connections and damage and reboot the system. Contact your dealer if the motor cable has suffered any visible damage.



Error number	Error description and solution proposal
3 & 4 <i>Error</i> Solution	Error Grip Offset <i>Throttle not in neutral position on system start.</i> Check that the throttle is in "neutral" position when you turn on your bicycle. Contact your dealer if the error does not disappear when the throttle is in neutral position – this means that it is probably defective.
5 <i>Error</i> Solution	Error Torque Offset <i>Torque sensor not installed properly or affected on system start.</i> Place your feet on the ground and reboot the system. If the fault does not disappear when there is no load on the pedals, the torque sensor is likely to be defective. Contact your dealer in order to get a replacement.
6 <i>Error</i> Solution	Error I Offset <i>Error in battery sensor, battery switch/connection is probably damaged.</i> Remove the battery and check the battery switches and connections. If the switches are in order, fit back the battery and reboot the system. Contact your dealer if any switches or connections are damaged.
7 <i>Error</i> Solution	Error Fast Over Voltage <i>Temporary overvoltage because of too high a speed or damaged battery switches or connections.</i> Remove the battery and check the battery switches and connections. If a switch or a connection is defective or damaged, contact your dealer for a replacement.
8 <i>Error</i> Solution	Error Slow Over Voltage <i>Lasting overvoltage.</i> Reboot the electrical system of the bicycle.
9 <i>Error</i> Solution	Error Fast Under Voltage <i>Lasting undervoltage. The battery is probably too weak or there are battery switches or connections that are damaged.</i> Remove the battery and check the battery switches and connections. If the switches are in order, fit back the battery and reboot the system. Contact your dealer if any switches or connections are damaged.
10 <i>Error</i> Solution	Error Slow Under Voltage <i>Lasting undervoltage. Battery is probably dead.</i> Recharge the battery and reboot the electrical system of the bicycle.
11 <i>Error</i> Solution	Error Over Temp Motor <i>Motor overheated.</i> Let the motor cool off and reboot the electrical system of the bicycle. Contact your dealer if the error does not disappear after 20 minutes.
12 <i>Error</i> Solution	Error Over Temp Controller <i>Motor controller overheated.</i> Let the motor controller cool off and reboot the electrical system of the bicycle. Contact your dealer if the error does not disappear after 20 minutes.



Error number	Error description and solution proposal
13 <i>Error</i> Solution	Error Parameter <i>General error.</i> Reboot the electrical system of the bicycle.
14 <i>Error</i> Solution	Error Under Temp Motor <i>Temperature too low for motor operation.</i> Bring the bicycle to room temperature – approx. 20°C. Contact your dealer if the error does not disappear after 1 hour.
15 <i>Error</i> Solution	Error EEPROM <i>Memory error.</i> Reboot the electrical system of the bicycle.
16 <i>Error</i> Solution	Error Parameter Property <i>Error in program management.</i> Reboot the electrical system of the bicycle.
17 <i>Error</i> Solution	Error Akku Temp <i>Battery is either too hot or too cold.</i> If overheated, let the battery cool off, and if too cold, store the battery at room temperature. Reboot the electrical system of the bicycle. Contact your dealer if the error does not disappear after approx. 1 hour.
18 <i>Error</i> Solution	Error Undefined Bike Constellation <i>System error.</i> Reboot the electrical system of the bicycle.
19 <i>Error</i> Solution	Error BMS <i>Battery error.</i> Remove the battery and recharge it. Reboot the electrical system of the bicycle.
20 <i>Error</i> Solution	Error BMS Version <i>Battery error.</i> Reboot the electrical system of the bicycle.
21 <i>Error</i> Solution	Error Torque Signal <i>Torque sensor not installed properly, cables/connections damaged or affected on system start.</i> Place your feet on the ground and reboot the system. Check the cables and connections for damage, in particular, around the pedals. Contact your dealer in case of visible damage
22 23 24 25 26 27 28 29 30 <i>Error</i> Solution	Error BMS Faulty high power current offset Error BMS Overcurrent discharge Error BMS BQ Short circuit Error BMS BQ Overload Error BMS BQ Overvoltage Error BMS BQ Undervoltage Error BMS BQ Overtemp Error BMS BQ Failed config Error BMS Faulty low current offset <i>Battery error.</i> Reboot the electrical system of the bicycle.



Error number	Error description and solution proposal
50 <i>Error</i> Solution	Error Display Communication <i>Communication error between controller and display.</i> Check that the connector between the display and the cable is installed properly and that the cable, connector or display is not damaged. Reboot the electrical system of the bicycle. Contact your dealer if there is any damage to the connector, display or cable – or if the error does not disappear.



13 EC DECLARATION OF CONFORMITY

Issued in accordance with EC Directive 93/42/EEC, as amended by Directive 2007/47/EC.

Manufacturer: PF mobility aps, Bjerregårdvej 15, Timring, DK-7480 Vildbjerg
Tel. +45 99 92 06 00
CVR No. 25059026

Models:

- Disco Junior
- Disco Senior

Versions: Electrically powered (250W) for all models:

- Direct power 36V

Product description: 3-wheeled electrically assisted bicycles and 3-wheel pedal-operated bicycles for users with special needs and handicaps.

Directives used: Directive 93/42/EEC concerning medical devices
and
Directive 2004/108/EC relating to electromagnetic compatibility
Class I

Equipment classification:
Harmonised standards used: EN 12182 – Assistive products for persons with disability –
General requirements and test methods.
EN 15194 – Electrically power assisted cycles
EN 61000-3-2, EN 61000-3-3, EN 61000-4-3, EN 61000-4-8
Electromagnetic compatibility.
EN 55014-1 + -2: Electromagnetic compatibility

Date of issue: 20.11.2015

PF mobility aps hereby declares that the products above meet all relevant provisions of Council Directive 93/43/EEC. The product is safe to use under conditions that can reasonably be expected, as described in the instruction manual.

This declaration is considered void for products that have undergone any kind of alteration without a prior written agreement with PF mobility aps.

PF mobility aps is obliged to set up procedures for continuous market monitoring and risk management in order to ensure corrective actions in the event of undesired events and risks.

Technical documentation for the products is kept by PF mobility aps at the address above.

Timring, 20.11.2015

Martin Søndergaard
CEO

