


## How to choose the correct size of a Theraplay foot propelled cycle.

In contrast to many children's bicycles, Theraplay cycles are sized by the leg length rather than the age as this is the most reliable way of getting the correct size of bike for each rider. Age or overall height alone may act as a guide but they are not a reliable measure as there is a large variance in average leg lengths and heights for each age. If you know the height of the rider, the following table will give an indication of the appropriate size of Theraplay cycle.


| Overall Height Of Rider:    | 70cm                           | 80cm | 90cm | 100cm | 110cm | 120cm | 130cm | 140cm | 150cm | 160cm | 170cm | 180cm | 190cm | 200cm |
|-----------------------------|--------------------------------|------|------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| IMP                         | 85cm (31/2") ----- 105cm (41") |      |      |       |       |       |       |       |       |       |       |       |       |       |
| Terrier/Terrier Hitch       | 100cm (39") ----- 125cm (49")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| TMX/TMX Hitch/TMX T5        | 120cm (47") ----- 145cm (57")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| Quad Small                  | 120cm (47") ----- 145cm (57")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| Tracer Jnr                  | 120cm (47") ----- 152cm (60")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| Tracker 20"/ Tracker T5 20" | 140cm (55") ----- 165cm (65")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| QUAD - Large                | 140cm (55") ----- 165cm (65")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| Tracker 24"/Tracker T5 24"  | 150cm (59") ----- 175cm (69")  |      |      |       |       |       |       |       |       |       |       |       |       |       |
| Tracer                      | 150cm (59") ----- 180cm (71")  |      |      |       |       |       |       |       |       |       |       |       |       |       |

For a more accurate guide to the finding the correct size of cycle take time to measure the rider's inside leg, this can be done in two ways:



1. Stand the rider in their socks against a wall. Place a small book between their legs, as high as is comfortable. Ask the rider to step away from the wall and measure from the top of the book to the floor.
2. The same measurement can be taken from a lying position. With the rider lying on their back, measure from the underside of the foot to the top of the inseam of the trouser leg. If the rider is not able to fully straighten their leg, measure to the maximum comfortable point of leg extension.

Once you have measured the inside leg, use the table below to find the right size of cycle for the rider. The seat height can be adjusted as the rider grows and the handlebar position can be adjusted to accommodate different arm reaches. PLEASE REMEMBER, THERAPLAY OFFERS AN ASSESSMENT SERVICE TO HELP YOU CHOOSE THE CORRECT MODEL TO SUIT THE INDIVIDUAL NEEDS OF EACH RIDER.



|                            | Overall Length   | Overall Width   | To Suit Inside Leg              | Seat to Handle Bars<br>Min - Max | Max User Weight |
|----------------------------|------------------|-----------------|---------------------------------|----------------------------------|-----------------|
| Imp                        | 1000mm (39 1/2") | 610mm (24")     | 320mm (12 1/2") - 430mm (17")   | 280mm(11") - 370mm(14 1/2")      | 32kg (70lbs)    |
| Terrier/Terrier Hitch      | 1140mm (45")     | 735mm (29")     | 420mm(16 1/2") - 545mm(21 1/2") | 360mm(14") - 430mm(17")          | 45kg (100lbs)   |
| TMX/TMX Hitch/TMX T5       | 1500mm (59")     | 735mm (29")     | 510mm(20") - 640mm(25")         | 480mm(19") - 560mm(22")          | 70kg (150lbs)   |
| Quad - Small               | 1450mm (57")     | 800mm (31 1/2") | 510mm(20") - 640mm(25")         | 480mm(19") - 560mm(22")          | 100kg (220lbs)  |
| Tracer Jnr                 | 1780mm (70")     | 800mm (31 1/2") | 510mm(20") - 710mm(28")         | 480mm(19") - 560mm(22")          | 70kg (150lbs)   |
| Tracker 20"/Tracker T5 20" | 1650mm (65")     | 800mm (31 1/2") | 610mm(24") - 740mm(29")         | 535mm(21") - 610mm(24")          | 90kg (200lbs)   |
| Quad - Large               | 1750mm (69")     | 800mm (31 1/2") | 640mm(25") - 760mm(30")         | 535mm(21") - 610mm(24")          | 100kg (220lbs)  |
| Tracker 24"/Tracker T5 24" | 1753mm (69")     | 800mm (31 1/2") | 640mm(25") - 760mm(30")         | 535mm(21") - 610mm(24")          | 90kg (200lbs)   |
| Tracer                     | 1930mm (76")     | 800mm (31 1/2") | 640mm(25") - 815mm(32")         | 560mm(22") - 640mm(25")          | 100kg (220lbs)  |